

INFORMATION TECHNOLOGY SERVICES 630.907.5000 . 1500 SULLIVAN ROAD, AURORA, IL 60506-1000 . IMSA.EDU

IMSA Verner Logger Pro Client

for Windows OS

TABLE OF CONTENTS

IMSA Verner Logger Pro Client	
TABLE OF CONTENTS	1
Downloading the Installation File(s)	1
Extract the Installation Files	2
Installing IMSA Verner Logger Pro Client	2

Downloading the Installation File(s)

- Go to https://imsa.edu/scs
- Click on the download link for the Installation/Setup for Windows OS
 - Science Class Requirement: Logger Pro
 - <u>Macintosh OS</u>
 - Windows OS
 - Instructions
 - <u>Download</u>
- Click the Dropbox download button

Download ~

• Make sure Save File is selected (if prompted) and click OK

,	What should Fire	fox do with this file?		
	Open with	Windows Explorer (default)	\sim
1	Save File			
h	Do this <u>a</u> uto	omatically for files like this fi	rom now on.	
			OK	Cancel

- Click Save to download the file to your Downloads folder on your Windows computer
- Close the pop-up windows and the web browser



INFORMATION TECHNOLOGY SERVICES

630.907.5000 . 1500 SULLIVAN ROAD, AURORA, IL 60506-1000 . IMSA.EDU

Extract the Installation Files

- Go to your downloads folder and open it so you can see the contents
- Find and right-click on the downloaded file (Install-LoggerPro.zip) and select Extract All from the pop-up menu

ivame	Date modified Type	
📕 Install-LoggerPro.zip	2/15/2022 12:56 DM Cpm	
	Open	
	Open in new window	
	Mount as ImDisk Virtual Disk	
	Extract All	
	Pin to Start	

 Click Extract to accept the default extract location. Note the location, as you will need it in the next step.



Installing IMSA Verner Logger Pro Client

- Navigate to the location noted above and double-click on the directory/folder to open it.
 - A window may automatically open to this location, depending on your Windows OS configuration.
- Double-click on the Install-LoggerPro folder
- Launch the installer, **installApplication (Run as Administrator).cmd**. Right-click and select Run as Administrator.

0 installApplication (Run as Administr	at> E/11/2020 2.20 DMA Mindaum Camana
	Onen
15 Logger Pro 3	open
LoggerPro3 12 Undate	Edit
Coggen too_nc_opdate	Duint
StopperPro3_14_Update	Print
I D202 and earlier Uninstall	Run as administrator
So croos_and_eanier_onnistan	· · · · · · · · · · · · · · · · · · ·
Les M	

- Click Yes when prompted to confirm running the installer.
- You will see a Windows protected your PC prompt. Click on the More Info Link



• Click the Run anyway button

Run anyway	Don't run

• You will see a User Account Control prompt, select Yes to allow the configuration script to run



INFORMATION TECHNOLOGY SERVICES 630.907.5000 . 1500 SULLIVAN ROAD, AURORA, IL 60506-1000 . IMSA.EDU

• The installer will run in the back ground for a few seconds with minimum activity on screen, you may see the following while installing the base package



• Once the base page has been installed, the update will be applied, you may see the following during the update process to version 3.16.2



- When the upgrade has complete the pop-up windows will disappear from your screen
- Installation Completed